

VICTUS88 COLLECTION GUIDE



Step 1:

Wash hands thoroughly with soap and water. Activate heat pack by bending the coin in the heat pack back and forth. You will hear a click and the gel inside will slightly change color.



Step 2:

Place heat pack vertically on upper arm, just below the shoulder, for 3 minutes. If heat pack does not warm up, rub the location with hand. The heat pack may be taken on and off if the temperature becomes uncomfortable.



Step 3:

Use the alcohol wipe to clean the heated area of your arm and let dry. Ensure entire heated area is wiped and allowed to dry for best results.



Step 4:

While your arm is drying, remove the collection device from its bag. Unscrew the tube cap and press the tube into the end of the device. Use a slight twist for a secure fit.



Step 5:

Remove the paper from the bottom of the device by pulling on the tab. This will expose the sticky side of the collection device.



Step 6:

Stick the device in the center of the heated portion of your arm with the tube facing forward. Failing to place the device in the center of the heated area may cause poor results.



Step 7:

Remove the safety clip by pulling it away from your body. Discard the safety clip.



Step 8:

Slide the lever on the device down until the audible “CLICK” is heard. Relax your arm by your side - there is no need to hold the device on your arm.



Step 9:

Fill tube per tube instructions, then remove the device by peeling it away from your arm. Remove the tube from the device using a gentle pull and twist.



Step 10:

Press the cap onto the tube firmly until a click is heard. A slight twist may be needed. Label the tube with your full name and date of birth. Wash hands after use.

